Instructions for Use: Endotracheal Tube

Not Made with Natural Rubber Latex

EN Endotracheal Tube

Federal (USA) Law restricts this device to sale by or on the order of a physician or any other practitioner licensed by the law of the State in which he/she practices to use or order the use of the device.

Intended Use: Endotracheal Tubes (cuffed or un-cuffed) are indicated airway management.

Description:

- Endotracheal Tube (cuffed or un-cuffed)
- Smooth Murphy Eye to provide an alternate ventilation port if tube becomes occluded.
- Softer rounded Murphy Eye is less invasive.
- Smooth tip to minimize trauma during intubation.
- Radiopaque strip embedded to assess exact location of tube.
- High volume, low pressure cuff provides an effective low pressure seal.
- Standard 15mm connector supplied with each tube, which insures compatibility with circuit connectors.
- Bold, easy-to-read depth marks for placement in low light conditions.
- Valve with sensitive pilot balloon indicates cuff's status.

Material: Non-irritating, Non-sensitizing, Non-cytotoxic

Contraindications: As with any PVC Endotracheal Tubes, use of Endotracheal Tubes in procedures which will involve the use of a laser beam or electrosurgical active electrode in the immediate area of the device is contraindicated.

Warnings/Cautions:

1. Verify proper assembly/attachment.
2. Test cuff by inflating it prior to use.
3. Deflate cuff prior to insertion.
4. Tube placement and intra-cuff pressure should be continuously monitored.
5. A second method of monitoring tube placement is recommended.
6. Safe cuff pressure should be less than or equal to 25 cm H₂O.
7. Reassess tube placement whenever patient movement, clinical interventions, or clinical condition suggest the possibility of tracheal tube dislodgement.
8. Beware of breathing circuit disengagement when transferring patients by stretchers.
9. Use proper tube fixation methods.
10. Do not use lubricants with a petrochemical base.
11. Single use only.
12. (Do not re-sterilize.) Resterilization could change the physical characteristics of the material and should not be attempted.
13. Do not use if package has been previously opened or damaged.
14. Exposure to elevated temperatures and ultraviolet light should be avoided during storage.
15. Store at a temperature that does not exceed 49°C.
16. Sterile if package is not opened or damaged.

Suggested Directions for Use:

1. Remove the sterile device from the package and retain package for patient record.
2. Test the cuff inflation system prior to intubation (in the case of cuffed Endotracheal Tube).
3. Deflate cuff prior to intubation (in the case of cuffed Endotracheal Tube).
4. If using a stylette, verify that the end of the stylette does not protrude from the end of the tube or pass through the Murphy Eye.
5. Inflate the cuff slowly while using the minimum amount of air required to provide an effective seal (in the case of cuffed Endotracheal Tube).
6. Monitor the cuff pressure levels during intubation in order to insure the effective seal (in the case of cuffed Endotracheal Tube).
7. Tracheal tube placement should be confirmed by viewing the position of the tube tip with chest radiograph.
8. Remove the tube only after complete deflation of the cuff (in the case of cuffed Endotracheal Tube).
9. Dispose device and device packaging in accordance to hospital practices and applicable state, federal and other regulations.

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