



Laryngoscope Recycling Program

Thinking beyond use.

When you use disposable laryngoscopes from SunMed, you are promoting better patient care by reducing the risk of cross-contamination. Additionally, you are eliminating the need for harmful chemicals that would normally be used for laryngoscope reprocessing. Thinking beyond use, SunMed supports the care of our environment through SunCycle, an easy and effective program to recycle disposable laryngoscope blades and handles. SunCycle is available to facilities who use SunOne and IntuBrite single use laryngoscope systems.

With SunCycle, you can help reduce waste, create a healthier community and reduce reprocessing costs for your facility.



Easy.

SunCycle Laryngoscope Recycling is simple:

- 1 Order SunCycle containers for recycling of any SunMed single use laryngoscope.
- 2 Following device use, dispose of designated SunMed products in the SunCycle receptacle.
- 3 When it is time for a new receptacle, call or email SunMed for a replacement.
- 4 When the new receptacle arrives, replace the full bin with the new empty one.
- 5 Package up the full receptacle and ship it back in the pre-labeled carton.



Category 'A' approved containers

The SunCycle Laryngoscope Recycling Program uses fully-compliant, Category A-approved, disposable waste containers. Available in a convenient 12-gallon size.

ROUND PORT HOLE
Easy disposal of laryngoscopes

HINGED LID SYSTEM
Not available elsewhere, and preferred by clinicians. Enables easier disposal of large volumes while keeping lid connected to the container.

16 POINT POSITIVE LOCKING LID SYSTEM
Audible clicks confirm the lid is locked.

PUNCTURE-PROOF, LEAK-PROOF LID
Competition is only puncture and leak-resistant. SunCycle containers feature leak-proof seals that permanently adhere the lid to the container –100% better than the traditional rubber gaskets on other PG-II containers.

INTEGRATED SIDE HANDLES
Makes it easy to move, even when full.

PROGRAM IDENTIFICATION
Clearly labeled with SunCycle recycling guideline.

TAMPER PROOF
Competition is only tamper resistant, but with SunCycle, no one can tamper with the contents of the container unless the top is cut off.

ENVIRONMENTALLY FRIENDLY
The containers burn clean; no toxic or harmful gas emissions.



SunCycle is intended for recycling of SunMed single use laryngoscope blades and handles. Other medical waste cannot be accepted. Do not place sharps, hazardous or radioactive-contaminated waste in this container system.



COMPONENT	ITEM	MAX CAPACITY	L x W x H
12 Gallon Bin	5-1080-12	55 lbs / 25 kg	13" x 17.3" x 20.8"

Your contribution.

When the laryngoscopes leave your facility, they arrive at a recycling center. There, the products are sterilized, sorted, disassembled, and processed to then re-enter the world as new products.

Compliant mailback program.

The SunCycle Program contains everything needed to properly and safely package and dispose of your SunMed disposable laryngoscope components. When properly packaged by the waste generator, this mail-back system meets all FedEx requirements for sending regulated medical waste through the U.S. mail to our disposal site. *Do NOT place this system in FedEx pick up if it is not properly and securely assembled.* Please read the Use, Shipping, and Packaging Instructions prior to using SunCycle in your facility.



SunCycle is a SunMed program
in partnership with



FOR MORE INFORMATION, CONTACT SUNMED
1-800-433-2797 | info@SunOneSystem.com
www.SunOneSystem.com/go-green/