



PRESSURE INFUSION BAGS

- Pressure relief valve prevents over inflation (330 mm Hg pressure relief)
- Large, oval shaped bulb allows for quick and easy inflation of bladder
- Dedicated fluid bag hook enables efficient loading and unloading of the system without having to remove the infusion bag from the IV pole
- Suitable for use with external inflation sources
- Color-coded gauge makes for accurate pressure monitoring (0-300 mm Hg)
- Incredibly reliable – 100% tested



Clear solution bag panel enables patient label and barcode scanning without having to remove bag.



Three-way stopcock ensures precise control of pressure while the large oval bulb allows for quick bladder inflation, all with ONE hand.



360° view of the color-coded gauge facilitates accurate pressure monitoring while the safety pressure relief valve prevents over-inflation.

ITEM	DESCRIPTION	SIZE
4005H	Infu-Surg	500 mL
4010H	Infu-Surg	1000 mL
4030	Infu-Surg	3000 mL
4050	Infu-Surg Clarity	500 mL
4100	Infu-Surg Clarity	1000 mL

PRODUCT DETAILS AVAILABLE AT WWW.SUN-MED.COM
FOR MORE INFORMATION CALL 1-800-433-2797