

ADJUSTABLE PEEP VALVE

LATEX FREE | SINGLE USE


- Designed for use with manual resuscitator/ventilator
- Pressure range: 5-20 cm H₂O
- Large print on scale for improved visibility of setting
- Connects to exhalation port of manual resuscitator




ITEM	DESCRIPTION	SIZE	PK
VP700	Adjustable PEEP Valve, 5-20 cm H ₂ O	19 mm	20
VP702	Adjustable PEEP Valve + Adapter, Converts to 30 mm	19 mm + 30 mm Adapter	20
VP705	MRI Conditional Adjustable PEEP Valve	19 mm	20

PEEP VALVE

LATEX FREE | REUSABLE



8-6002-30



8-6003-30

- 5-20 cm H₂O
- Designed for use with manual resuscitator/ventilator
- Connects to exhalation port of manual resuscitator

ITEM	DESCRIPTION	SIZE	PK
8-6002-30	5-20 cm H ₂ O PEEP Valve	30 mm ID x 33 mm OD	1
8-6003-30	PEEP Valve Diverter	30 mm ID	1