

GREENLINE®/D™ FIBER OPTIC ALL-METAL BLADE


LATEX FREE | STERILE | SINGLE PATIENT USE | DISPOSABLE | CE MARK

Fiber optic, all-metal laryngoscope blades solve contamination problems and eliminate the cost and time spent cleaning blades and returning them for use. These surgical grade 303/304 stainless steel disposables also answer the professional's request for a non-plastic disposable blade. The flexible fiber optic bundle is protected in a black plastic sheath. All sizes are compatible with any reusable "Green System" Fiber Optic handle and comply with ISO 7376 Green Standard.



ALL-METAL MILLER LARYNGOSCOPE BLADES

ITEM	DESCRIPTION	SIZE	A	B	PK
5-5336-20	Premature	00	67 mm	10 mm	20
5-5336-00	Neonate	0	77 mm	10 mm	20
5-5336-01	Infant	1	92 mm	10 mm	20
5-5336-15	Small Child	1.5	130 mm	12 mm	20
5-5336-02	Child	2	100 mm	12 mm	20
5-5336-03	Medium Adult	3	130 mm	14 mm	20
5-5336-04	Large Adult	4	155 mm	18 mm	20



ALL-METAL MACINTOSH LARYNGOSCOPE BLADES

ITEM	DESCRIPTION	SIZE	A	B	PK
5-5335-00	Neonate	0	80 mm	16 mm	20
5-5335-01	Infant	1	92 mm	22 mm	20
5-5335-02	Child	2	100 mm	21 mm	20
5-5335-03	Medium Adult	3	130 mm	22 mm	20
5-5335-35	Ext. Med. Adult	3.5	145 mm	25 mm	20
5-5335-04	Large Adult	4	155 mm	25 mm	20

