

CapnoVue® Scope Convenience Kits

LATEX FREE | SINGLE PATIENT USE | DISPOSABLE | DEHP FREE

Simplify life in your busy endoscopy lab with the Salter Lab® Blue Series Kits. The CapnoVue® Scope Convenience Kit combines our scope-ready mask with your choice of a standard Luer-Lok® or reflective connector EtCO₂ sample line.




47MSF8-7-6-50
CapnoVue Scope Mask
plus EtCO₂ Sample Line



CAPNOVUE® SCOPE CONVENIENCE KITS

- CapnoVue masks provide higher O₂ flow rates with oral-nasal sampling
- Flow dynamics separate gasses by pulling CO₂ toward vents and sampling lines and away from O₂ port
- Each available in two tubing lengths



ITEM#	DESCRIPTION	MASK O ₂ LENGTH	ETCO ₂ LENGTH	CONNECTOR	CS
47MSF8-7-6-50	CapnoVue Scope Convenience Kit	7' (2.1 m)	6' (1.8 m)	Reflective	50
47MSF8-14-13-50	CapnoVue Scope Convenience Kit	14' (4.3 m)	13' (4 m)	Reflective	50
4745-7-7-50	CapnoVue Scope Convenience Kit	7' (2.1 m)	7' (2.1 m)	Luer-Lok	50
4745-14-15-50	CapnoVue Scope Convenience Kit	14' (4.3 m)	15' (4.6 m)	Luer-Lok	50

CapnoVue is a registered trademark of Salter Labs. Luer-Lok and Luer-Slip are registered trademarks of Becton Dickinson.