


CO₂ EASY® CARBON DIOXIDE DETECTOR

LATEX FREE | SINGLE PATIENT USE | DISPOSABLE | DEHP FREE

With a large two inch visual indicator, the CO₂ Easy® makes it easy to view critical carbon dioxide levels without looking away from the patient. The nontoxic, pH sensitive indicator uses familiar purple-to-yellow color conversion when CO₂ is detected.



- Designed for fast, easy airway assessment.
- Verifies endotracheal tube placement.
- 3-year shelf life.
- Complies with Standard of Care Guidelines from the American Heart Association, American College of Emergency Physicians, American Association for Respiratory Care and the American Society of Anesthesiologists.

CO₂ EASY® CARBON DIOXIDE DETECTOR

ITEM	DESCRIPTION	SIZE	WEIGHT	PK
562134	CO ₂ Easy	Adult	> 15 kg (33 lbs)	10
562135	Pedi CO ₂ Easy	Pediatric	1 - 15 kg (2.2 - 33 lbs)	10